


FloorBed[®] 2





FloorBed[®] 2
By Accora

FloorBed[®] 2

Reduce
Injury

Reduce
Carer
Burden



Manage
Risk



TECHNICAL SPECIFICATION

Dimensions:	94 cm x 228 cm (37" x 90")	Lifting range:	7 cm - 80 cm (2.75" - 31.5") <i>(excl. mattress)</i>
Platform size:	90 cm x 200 cm (36" x 79")	Profiling functions:	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Safe working load:	185 kg (29 st)	Part number:	NSB-0-FL2-200
Max patient weight:	150 kg (24 st)		

How does FloorBed® 2 compare to a traditional nursing bed?



FLOORBED® 2 FEATURES AND FUNCTIONS:

FloorBed® 2 raises to 800 mm which facilitates an even safer working height for taller carers

The height of the head and foot board when FloorBed® 2 is at full height is 110 cm

The low head and foot board enable access to care tasks such as hair care and lower limb dressing

Low head and foot board will facilitate access to sliding sheets at the head or foot of the bed allowing single carer use to reposition a person up or down the bed

FloorBed® 2 Features



Width adjustable mattress guides



No trailing wires



Full nursing height



Concealed wheels



Backrest and kneebreak functions



Mattress platform height of just 7.1 cm – similar to a tennis ball!

FloorBed® 2 Accessories



Integrated side rails made of soft material



Head and foot board bumpers reduce impact and minimise the risk of injury*



Safety mat, 22 cm high. Provides same height as a FloorBed® with 15 cm foam mattress*



Safety sleeve reducing risk of entrapment as the platform mechanism is protected*



Bed lever to assist with bed mobility



Rotating bed lever for extra support in movement



Support bar to assist with bed mobility




Easy to transport

* The picture depicts FloorBed® 1


Visit our website at www.floorbed.com for the full range of FloorBed® 2 accessories.

Why FloorBed®?

NURSING AND LOW BEDS



FLOORBED®



When to use FloorBed®?

If the user:

- Is at risk of bed falls
- Has uncontrolled movements
- Has difficulties with leg lifting, difficulty with bed transfers and difficulty with bed mobility



Injury risks when falling from bed

